VENTA-NEB[®]-ir A-I-C-I[®]


Microprocessor Controlled • Mobile Ultrasonic Nebulizer
for VENTAVIS[®] Inhalation

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Operating Instructions

für Heimtherapie/spontan atmende Patienten
for home therapy / spontaneously breathing patients

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VENTA-NEB®-ir A-I-C-I®

Microprocessor Controlled - Mobile Ultrasonic Nebulizer

Operating Instructions

Unit Type VN-100/4-2.4 MHz

Made in Germany

Dear Patient,
With the mobile Ultrasonic Nebulizer VENTA-NEB®-ir you received an inhalation device that was adjusted extremely precise in a conditioned room.
In order to ensure a constant operation of the equipment we ask you to carefully read the operating manual and to follow the instructions before you put the device into service.
We wish you every success for your treatment with the VENTA-NEB®-ir.

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Stempel/Stamp:

Table of Contents

1.0	IEC Symbols
2.0	Important Safeguards
3.0	Intended Purpose
3.1	Explanations on Operating Modes/Nebulization Programs
3.2	Functional Description - Spontaneous Respiration/Home Therapy
3.3	Aerosol Spectrum (Particle Sizes)
4.0	Important Parts of Your Ultrasonic Nebulizer
5.0	How to Operate Your Ultrasonic Nebulizer Using the Mouthpiece
6.0	Power Supply of the Ultrasonic Nebulizer
7.0	Using Your Ultrasonic Nebulizer
7.1	Programming and Settings
7.1.1	VENTA-NEB®-ir 2,4 MHz
8.0	Using an Extension Hose for Inhalation when Lying
8.1	Using a Face Mask for the Inhalation (for Children)
9.0	Cleaning
9.1	Autoclaveable Plastic Parts
9.2	Contact Fluid Chamber and Control Unit
10.0	Changing the Medicine Cup
10.1	Replacement Intervals for Filter Membrane/Medicine Cup
11.0	Maintenance
12.0	Notes on Troubleshooting
13.0	Specifications
14.0	Accessories / Order Information VENTA-NEB®-ir
15.0	Warranty
16.0	Declaration of Conformity
17.0	Warranty Card (tear-off card)

1.0 IEC Symbols

Attention, consult operating instructions

Protection class II equipment

Type B application part

Class AP equipment

1275

2.0 Important Safeguards

Each handling of the device requires the exact understanding and observance of the present operating instructions.

In any case, the liability for the safe functioning of the device passes over to the user, if intervention by third persons or handling not according to the intended use occur.

Important information is highlighted by these terms:

DANGER

Urgent safety information for hazards that might cause serious injury.

WARNING

Important information for operating steps that might cause malfunctions of the device.

CAUTION

Information preventing damage to the product.

NOTE

Information to which you should pay particular attention.

[Components – see Sketch 4.0]

Please read the operating instructions carefully before the first operation of the product. Keep the operating instructions in a safe place.

DANGER

1. Always unplug the device immediately after use.
2. Do not use the device while taking a bath.
3. Do not place the device where it can fall into water.
4. Do not immerse the device in water or other liquids.
5. Do not use the device when it has fallen into water. Unplug it immediately.
6. Do not use the device under rain.
7. Do not use the device near flammable agents
8. Never insert hands or fingers into the contact fluid chamber or the medicine cup while the device is operating.

WARNING

1. Activated HF-Communication equipment (like cellular phones, etc.) nearby the **VENTA-NEB[®]-ir** may affect its proper function.
2. Electric equipment should never be left unattended when plugged in.
3. Close supervision is necessary when the device is used by or near children or invalids.
4. Only use this device for the intended purposes as described in this operating instructions.
Do not use accessories not recommended by the manufacturer.
5. Never operate this product if
 - a) the power cord or the plug are damaged.
 - b) the device is not working properly.
 - c) the device has been dropped or damaged.
 - d) the device fell into water. In such cases return the device to the manufacturer or an authorised NEBU-TEC service dealer for inspection and repair.
6. Keep the power cord away from heated surfaces.
7. Place the device on a flat and solid surface in such way that no air openings are blocked.
8. Never use the device while you are sleeping.
9. Never insert any object into the openings of the equipment.
10. Never use the dish washer or microwave oven for cleaning the Ultrasonic Nebulizer (never expose the device to direct microwave radiation).

11. When cleaning the contact fluid chamber pay close attention to that no moisture can penetrate the housing from outside.
12. Change the contact fluid (distilled water) every 24 hours.
13. Do not use any cleaning solutions, vinegar, hot or even boiling water, etc... for cleaning the **control unit, the contact fluid chamber, the transducer or the sensor**.
14. Do not use cleaning solutions on aqueous alkaline basis, aromatic hydrocarbons, ammonia and amines for cleaning the **autoclaveable plastic parts**. Use cleaners on saturated aliphatic hydrocarbon basis, alcohols, dilute mineral acids, neutral and acid salt solution instead.

NOTE

While the **VENTA-NEB®-Ir** is in operation no measures need to be taken regarding the electromagnetic compatibility.

Never lend the device to third persons.

Only inhale the medication prescribed by the physician.

Upon longer use the device may heat up on the bottom side.

3.0 Intended Purpose

Your **VENTA-NEB®-Ir** Ultrasonic Nebulizer is a portable device designed to deliver aerosols of different particle sizes (see section 3.3). This ensures an optimum deposition of the medication.

3.1 Explanation on Operating Modes/Nebulization Programs

In order to allow the user to control which program is set and activated the active program briefly (for app. 1 second) appears in the display after the Ultrasonic Nebulizer was connected to the electric circuit.

- **VENTA-NEB®-Ir** 2,4 MHz: 2 Programs (P1 or P2)

3.2 Function of A-I-C-I® for spontaneously breathing patients/home therapy

A-I-C-I® (active intermitted controlled inhalation)

When the device has been put into operation it indicates when and how often the patient should breathe in (inhale).

Acoustic signal: Exhalation

Acoustic and optic signal (green light [A11.1]): Inhalation

Due to this inhalation-scheme, a more efficient and a precise dosage can be guaranteed.

The produced vapour can be inhaled through the mouthpiece (or mask). The exhalation also takes place over the mouthpiece [A2] or the face mask and is controlled by valves installed inside the filter shell [A8, A9] that cannot be interchanged. These filters prevent any leakage of aerosol into the air and thus form a closed system of the nebulization unit.

Located on top of the Dome [A4] of device is a Luer/Lock Connection [A6] where the medication can be added without opening the nebulization unit.

This connection can also be used to add oxygen.


3.3 Aerosol Spectrum (Particle Sizes)

The particle size (MMAD(Mass Median Aerodynamic Diameter in μm) of the aerosol using

the green baffle plate and VENTAVIS[®] solution is 2.3 μm ,

the blue baffle plate and VENTAVIS[®] solution is 3.2 μm .


Baffle Plate 1 (VN-117G)	green colour	MMAD 2.3 μm measured with VENTAVIS [®]
Baffle Plate 2 (VN-117B)	blue colour	MMAD 3.2 μm measured with VENTAVIS [®]

NOTE

The VENTAVIS-ir can be adjusted to patient-specific/individual inhalation patterns. However, only the manufacturer or the treating physician are authorized to perform such adjustments.

4.0 Important Parts of Your Ultrasonic Nebulizer VENTA-NEB®-ir


- A1** VENTA-NEB®-ir Ultrasonic Nebulizer
- A2** Mouthpiece (Art. VN-105)
- A3** Exhalation Piece (Art. VN-104)
- A4** Dome (Art. VN-103) with Sealing Ring (Art. VN-110), Baffle Plate (Art. VN-117B/G/R/Y) and Luer/Lock Protecting Cap (Art. VN-116)
- A5** Rest for Luer/Lock Protecting Cap
- A6** Luer/Lock Connection
- A7** Luer/Lock Protecting Cap (Art. VN-116)
- A8** Exhalation Filter Shell with Valve (Art. VN-101)
- A8.1** Exhalation Filter Membrane (Art. VN-109)
- A9** Inhalation Filter Shell with Valve (Art. VN-101)
- A9.1** Inhalation Filter Membrane (Art. VN-109)
- A10** Medicine Cup (VN-102)
- A11** Multifunctional Indicator Light
- A11.1** Light as a request to inhale
- A12** Display
- A13** Infrared Sensor
- A14** On/Off Button
- A15** Start/Stop Button
- A16** Rechargeable Battery Pack (Art. VN-MCA)
- A17** AC Power Adapter 110-230/230 VAC (Art. VN-100N)
- A18** 12 VDC Car Power Cord for Cigarette Lighter Socket (Art. VN-100Z)
- A20** Tube Extension 22 OD/15 ID (Art. VN-B-109)
- A21** Face Mask for Children with Expiration Valve, Sizes 1/2/3 (Art. VN-122/123/124)
- A22** Luer/Lock Tube for Oxygen or Medicine Nebulization Control Line (Art. VN-111)
- A26** Adapter for adaptation of face mask 22 OD/22 OD - 15 ID (Art. VN-119)
- A27** Inhalation trainer

5.0 How to Operate Your VENTA-NEB®-ir Ultrasonic Nebulizer (USN) upon Inhalation through the Mouthpiece

Preparing the Ultrasonic Nebulizer


- 5.01** Remove the power plug [A17] so that the nebulizer will not start running unintentionally.
- 5.02** Fill the contact fluid chamber up to the blue mark with distilled or demineralised water using the measuring cup. The sensor must be entirely covered with this contact fluid (app. 45 ml).

WARNING

The utilization of other contact fluids (e.g. tap water, sterile water or saline) is strictly prohibited as this may result in an essential affection of the device's performance and even the total failure of the device.

- 5.03** Insert the medicine cup [A10]. Please ensure that the tip of the medicine cup is submersed in the contact fluid.


NOTE

The medicine cup [A10] is a disposable article. Inspect it for any defect before each use.


If the medicine cup is damaged or the unit is experiencing low output, replace the medicine cup or clean the transducer.

- 5.04** Make sure that the baffle plate is properly fastened inside the dome [A4], a sealing ring is inserted, and the Luer/Lock protecting cap [A7] is connected to the dome.


- 5.05** Insert the filter shell with the inhalation filter membrane [A9] in the designated opening of the dome [A4].

- 5.06** Place the dome [A4] with the inhalation filter [A9] oriented towards the back of the unit, and turn it clockwise until you hear a clear click. This sound is generated when the medicine cup [A10] connects to the dome [A4].


- 5.07 Assemble the mouthpiece [A2], the exhalation piece [A3], the exhalation filter shell [A8] and connect them to the dome [A4].
- 5.08 Fill-in the fluid (medicine solution) to inhale through the designated opening [A6] (Luer/Lock connection) in the top of the dome [A4] by means of a syringe.
- 5.09 Close the Luer/Lock connection [A6] with the protecting cap [A7].
- 5.10 Put the device into operation as described in section 7.0 – 7.1 of this operating instructions.

ENGLISH


NOTE

Upon several inhalations per day or when the device is transported between the single inhalations, the dome [A4] should be closed with the delivered Luer/Lock protecting cap [A7] and the two plugs for hygienic reasons. Then the mouthpiece [A2], the exhalation piece [A3] and the two filter shells [A8, A9] have to be stowed in the safety box.

6.0 Power Supply of the Ultrasonic Nebulizer VENTA-NEB®-ir

You can connect the Ultrasonic Nebulizer to three different power sources: alternating current - 110/220 volts, direct current - 12 volts (car), or rechargeable battery pack.

6.1 AC Operation

Plug the AC Power Adapter [A17] in an AC power outlet and the other end in the unit (110 or 220/230 volts).

NOTE for AC Operation

Do not use the device while taking a bath.
Do not place the device where it can fall into water.
Do not immerse the device in water or other liquids.
Do not use the device when it has fallen into water.
Immediately unplug the power adapter.

6.2 DC Operation

Plug the car power cord [A18] in the cigarette lighter socket and the other end in the unit.

NOTE for DC operation

Do not use the device while taking a bath.
Do not place the device where it can fall into water.
Do not immerse the device in water or other liquids.
Do not use the device when it has fallen into water.
Immediately unplug the car power cord.

6.3 Instruction guide of the rechargeable battery pack

Please find enclosed the separate instructin guide.

6.3.1 For charging purposes please only use the original power cord [A17]. Therefore, please connect the plug of the power cord to the rechargeable battery pack.

6.3.2 Depending on the charging-state, the rechargeable battery pack takes up to 12 hours until it is completely charged.

6.3.3 Please do not use the rechargeable battery pack while it is being charged.

6.3.4 Shortly after the charging process has been fnished, the letter combination „E“, „n“, „d“ (= Ende) appears successively on the display. Now the rechargeable battery pack is ready for use.

6.3.5 If you do not want to use the rechargeable battery pack immediately after charging it, it can remain connected to the power cord for a period of time (even a few days) although the charging process has been finished.

- 6.3.6 Please only connect the rechargeable battery pack to the **VENTA-NEB®-ir** while executing the inhalation.
- 6.3.7 Please remove the plug of the rechargeable battery pack from the nebulizer after each inhalation.
- 6.3.8 In order to check the charging-state, please push the button (Akkutest) for approximately 3 seconds. Shortly after having pushed this button, the charging capacity is shown on the display in per cent (%).

ENGLISH

Please ignore the green, red or orange LED display on the battery pack itself, since these indicator lights do not have any significance concerning the actual state of charge of the battery pack.

When the rechargeable battery pack is completely charged, it is sufficient for approximately one week.

NOTE for battery operation

Do not use the device while taking a bath.
 Do not place the device where it can fall into water.
 Do not immerse the device in water or other liquids.
 Do not use the device when it has fallen into water.
 Immediately unplug the rechargeable battery pack.

CAUTION

To prevent damage to the Ultrasonic Nebulizer and to ensure compliance with the EMV EN 55011 guidelines, only the original AC power adapter [A17] or the original rechargeable battery pack [A16] may be used.

NOTE


Bring defective battery packs to battery disposal points for proper disposal, or send them back to NEBU-TEC GmbH.

7.0 Using Your Ultrasonic Nebulizer VENTA-NEB®-ir

After connecting the **VENTA-NEB®-ir** to the power circuit the inhalation program set by your physician briefly appears in the display [A12] (for app. 1 sec.).

7.01 In order to turn on the ultrasonic nebulizer, please push the button „On/Off“ – a yellow light appears on the nebulizer. The amount of inspiration cycles is now shown on the display.

7.02 In order to commence the inhalation, please push the button “start/stop” [A15] (accompanied by an acoustic signal) and

breathe out slowly. Please breathe in regularly and slowly when noticing the next acoustic signal which is accompanied by an optic signal at the same time [A11.1] (green light appears). After breathing in, please shortly stop breathing and then breathe out slowly. (In order to learn how to inhale correctly and to optimize the inhalation, please use the “inhalation coach” [A27]). Please connect the „inhalation coach“ to the inhalation filter shell.

7.03 To interrupt the aerosol production, please press the **Start/Stop** sensor button [A15]. The multifunctional indicator light [A11] is yellow and the display [A12] shows PA (Pause). To continue the aerosol production press the **Start/Stop** sensor button [A15] again.

The display will change from "PA" to the time mode, the multifunctional indicator light is green, the nebulization is continued.

The inhalation is terminated when En (End) appears in the display. At the same time an acoustic signal sounds at the end of the inhalation.

7.04 Place the mouthpiece [A2] in your mouth and inhale the medicated aerosol over the inhalation filter and the valve by taking a slow, deep breath. The exhalation also takes place over the mouthpiece and the exhalation filter with valve.

The inhalation instruction: ‘The Right Way to Inhale’ is available separately.

7.05 Continue to inhale until the respective nebulization program has expired. This is indicated by an acoustic signal and En (End) in the display. The inhalation period can depend on the respective nebulization program as well as the breathing technique of the patient.

7.06 Switch off the device after the inhalation is finished (**On/Off**) [A14] and disconnect the Ultrasonic Nebulizer from the respective power source.

NOTE

The device is equipped with a multifunctional indicator light which indicates the operating state.

ADVICE

The device has a multi-function light [A11] and a light as a request to inhale [A11.1].

Light as a request to inhale [A11.1]:

Green light is on - Inhalation starts.

Green light is off - end of the inhalation or pause.

Multifunctional indicator light [A11]:

Yellow light - Device ready for operation

Green light - Device in operation

Red light - Malfunction

Display readings [A12]

(LB) Low Battery...

(LH) Low Hydrogen - contact fluid missing, or wrong contact fluid in the chamber

(SA) Impure or saline fluid in the unit (tap water, sodium chloride, mineral water, etc.)

(PA) Pause

(En) End...

NOTE


As the proportion of active substance of your medicine within the remaining residue is too small, please do not use it again for a new inhalation.

In order to check the amount of residue you may turn the device upside down and read the residue at the graduation on the dome [A4].

The device is leak-proof up to 7,5 ml of liquid in the medicine cup.

NOTE

The residue of the medication remaining in the medicine cup has to be emptied after each inhalation. In order to empty the residue remove the inhalation filter [A9] and pour out the remaining medication by tilting the device.

7.1 Programming and Settings of the Ultrasonic Nebulizer

The selection of the nebulization programs has to be conducted by the attending physician and authorized technical personnel that works on the order of the physician's recommendation.

The manufacturer of the **VENTA-NEB[®]-ir** does not issue any recommendations regarding the dosage of medication. For the inhalative application please carefully read the enclosed leaflet of the medicine.

7.1.1 Programming and Settings of the VENTA-NEB[®]-ir 2.4 MHz

The **VENTA-NEB[®]-ir 2.4 MHz** Ultrasonic Nebulizer allows the patient to choose between 2 different nebulization programs.

After connecting the **VENTA-NEB[®]-ir 2.4 MHz** to the electric circuit the set nebulization program is briefly (app. 1 second) shown in the display.

P1 Program 1 5,0 µg active substance on the mouthpiece 25 inhal. cycles

P2 Program 2 2,5 µg active substance on the mouthpiece 10 inhal. cycles

Changing the Nebulization Program on the VENTA-NEB[®]-ir 2.4 MHz

To change the programs on the **VENTA-NEB[®]-ir 2.4 MHz** please proceed as follows: Press both display buttons (Start/Stop and On/Off) [A14, A15] at the same time, the program shown in the display [A12] flashes (P1 or P2).

Set the desired program (P1 or P2) by pressing the button Start/Stop [A15] (down) or the button On/Off [A14] (up). 10 seconds after this procedure the display [A12] stops flashing and the **VENTA-NEB[®]-ir 2.4 MHz** shows the selected program.

For the operation of the device refer to the description in Chapter 7.0.

Setting:

Pressing both sensor buttons at the same time:


Display blinks

Pressing left button **Start/Stop** [A15]:

Set lower value.

Pressing right button **On/Off** [A14]:

Set higher value.

ENGLISH

If no button is pressed the set value will be stored after app. 10 seconds.


NOTE

The desired setting of the device is activated by one (short) tap on the sensor buttons [A14, A15]. If the buttons are pressed or held for a longer time this results in the respective setting being activated/deactivated/activated/deactivated...(on/off/on/off...) and so on.

Therefore we would like to ask you to actuate the two sensor buttons by shortly tapping on them.

8.0 Using an Tube Extension for Inhalation when Lying


An extension hose or tube extension [A20] can be used if the patient is lying while using the device. Then the hose [A20] is inserted between the mouthpiece [A2] and the exhalation piece [A3].

8.1 Using a Face Mask for the Inhalation (for Children)

A facial mask [A21] can be used especially for children. The face mask can either be directly connected to the exhalation piece [A3] or to a tube extension [A20], if necessary.

Then the face mask is connected to the exhalation piece [A3] or to the extension hose [A20] by means of an adapter.

9.0 Cleaning

Carefully following the steps outlined below will help maximise the performance and extend the service life of your Ultrasonic Nebulizer **VENTA-NEB®-lr**.

WARNING

To prevent possible risk of infection from contaminated accessory parts, we recommend to follow the manufacturer's instructions.

DANGER

Please disconnect the device from any power source before cleaning your Ultrasonic Nebulizer.

9.1 Autoclaveable Plastic Parts

Cleaning the Dome [A4] and Baffle Plate(s), Exhalation Piece [A3], Filter Shells [A8/A9] and Mouthpiece [A2].

The parts described above have to be cleaned upon inhalation once a day after the inhalation, or upon change of medication or inhalation several times a day after the last inhalation.

The parts are temperature resistant up to 134°C and have to be cleaned after 24 hours at the latest as follows:

How you should clean your device at home:

9.11 Disassemble the Dome [A4] with Sealing Ring and Baffle Plate(s), Exhalation Piece [A3], both Filter Shells [A8/A9] and Mouthpiece [A2]. You can push out the baffle plate through the opening in the top of the dome.

- 9.12 Clean the accessories under warm tap water or in the dishwasher every day.
- 9.13 After cleaning, scald out the plastic parts in a pot for app. 10 min.
- 9.14 Then allow the component parts to air dry. (Towel drying could lead to contamination or soiling.)
- 9.15 If using the microwave vapour sterilizer, please follow the separate operating manual. The Luer/Lock protecting cap [A7] must **N O T** be cleaned or sterilized in the microwave oven.
- 9.16 Carefully reassemble all parts. If the baffle plate has been removed from the dome [A4], please check that it was properly reinstalled.
- 9.17 All transparent parts [A2, A3, A4, A7, A8] are autoclaveable up to 134°C.

WARNING

Never put the Ultrasonic Nebulizer **VENTA-NEB®-ir** [A1] in a microwave oven.


IMPORTANT NOTE

The Dome [A4] with Sealing Ring and Baffle Plate(s), Exhalation Piece [A3], both Filter Shells [A8/A9] and Mouthpiece [A2] as well as Luer/Lock protecting cap [A7] should be replaced after 3 months when the unit is used several times a day. If used only once a day the parts mentioned above have to be replaced depending on wear and hygienic condition.

9.2 Contact Fluid Chamber and Control Unit of the VENTA-NEB®-ir**DANGER**


Never submerge the control unit of the Ultrasonic Nebulizer in water or cleaning solution.

- 9.21 Before cleaning the control unit, disconnect the AC power adapter [A17] from the device.
- 9.22 Only wipe the exterior of the control unit with a damp cloth.
- 9.23 Empty the contact fluid chamber after the last daily inhalation. Then carefully dry the inside of the contact fluid chamber with a soft cloth. After the cleaning turn the device upside down (place it on an absorbent pad) and let it dry in this position until the next utilization.

9.24 The transducer (at the bottom of the contact fluid chamber) should be cleaned once a week by wiping carefully with a cotton swab (performing rotating movements).

Top view of the USN


Clean transducer with a cotton swab performing rotating movements

CAUTION

Do not scratch the transducer with sharp-edged items.
Never push too hard on the transducer (at the bottom of the contact fluid chamber). Doing so may lead to damages.

10.0 Changing the Medicine Cup [A10]

If you remove the dome [A4] from the device, the medicine cup [A10] is attached to the dome by four lugs. The medicine cup is removed from the dome by pressing on them laterally.

10.1 Replacement Intervals for Filter Membrane and Medicine Cup

10.10 In order to ensure a secure dosing the medicine cup has to be replaced every day.

10.11 The exhalation filter membrane [A8.1] has to be changed every day.


10.12 Upon inhalation several times a day it might be necessary to change the filter membrane [A8.1] more often (increase of the exhalation resistance of the filter membrane).

10.13 The inhalation filter [A9.1] has to be changed once a week.

10.14 Turn the filter shell [A8/A9] anticlockwise to open it.


10.15 Replace the filter membrane [A8/A9].

10.16 Turn the filter shell [A8/A9] clockwise to close it again.

NOTE

The medicine cups are disposables and have to be replaced every day for hygienic and technical reasons (see section 10.1). Failure to observe the provided replacement intervals may lead to deformations of the medicine cup [A10]). These micro deformations of the medicine cup [A10] may considerably reduce the output of the Ultrasonic Nebulizer

ENGLISH

11.0 Maintenance

The **VENTA-NEB®-ir** Ultrasonic Nebulizer should be serviced every 2 years. All maintenance must be performed by NEBU-TEC GmbH or an authorised and qualified NEBU-TEC service dealer.

DANGER

Do not remove the exterior cabinet of the control unit. Non-compliance will lead to loss of warranty.

12.0 Notes on Troubleshooting

If you think your device is not working properly, please take a few moments to check for and repair these possible causes before you complain about it.

Symptoms	Possible Causes	Remedies
Display  (LB/Low Battery)	1. AC adapter defective. 2. Battery empty.	1. Contact manufacturer or dealer 2. Charge battery pack.
Display  (LH/Low Hydrogen)	1. No contact fluid in the contact fluid chamber. 2. Sterile or too pure water in the contact fluid chamber	1. Fill in 45 ml of contact fluid (sensor must be covered). 2. Add app. 1 ml of tap water to the 45 ml of contact fluid.

Display

(SA/Salt Recognition)

1. Saline or impure fluid in the contact fluid chamber (e.g. tap water, NaCl, mineral water).

1. Carefully wash out the chamber several times with distilled water. Carefully clean the sensor in the contact fluid chamber with a cotton swab or something alike, rinse out the chamber again with distilled water, and then refill the contact fluid chamber.

Reduced aerosol output (too much residue)

1. Worn or damaged medicine cup.
2. Contact fluid level in the contact fluid chamber too high/low.
3. Contact fluid chamber was not cleaned properly.
4. Multiple medicine cups placed in the contact fluid chamber.

1. Replace medicine cup.
2. Fill contact fluid chamber with distilled water up to the mark (app. 45 ml).
3. Clean device according to instructions.
4. Only place one medicine cup.

Device does not produce aerosol

1. Multiple medicine cups placed in the contact fluid chamber.
2. Worn or damaged medicine cup placed in the chamber.
3. Device not connected to power source.
4. No contact fluid filled in the contact fluid chamber.
5. No fluid (medicine solution) in the medicine cup.

1. Insert only one medicine cup.
2. Insert new medicine cup.
3. Connect device to power source.
4. Fill contact fluid chamber to proper level.
5. Fill medicine cup.

Inhalation or exhalation is more difficult

1. Filter membrane is clogged or soaked.
2. Dome is not properly fastened.

1. Replace filter membrane.
2. Check if dome is properly fastened.

13.0 Specifications VENTA-NEB®-Ir

Size	98 x 66 x 105 mm
Weight, control unit	280 g
Types of power supply	110/230 VAC power adapter
	12 VDC car power cord (cigarette lighter)
	12 VDC battery pack
Power supply	12 VDC, 1.5 A maximum
Operating power consumption	18 Watt maximum
Ultrasonic frequency	2.4 MHz (nominal)
MMAD	2.3 µm (green coloured baffle plate)
Medicine cup capacity	7.5 ml maximum
Contact fluid chamber capacity	45 ml
Electric protection class	II Type B

14.0 Accessories / Order Information VENTA-NEB®-Ir

Article Number	Designation	Packing Unit/Quantity
VN-100/2	VENTA-NEB®-Ir 2.4 MHz	1
VN-MCA	Microprocessor controlled rechargeable battery pack ON-2000	1
VN-100Z	12 VDC Car Power Cord (Cigarette Lighter)	1
VN-100N	AC Power Adapter FW 7555M/12 110/230 VAC	
	(available with international Adapters upon request)	1
Non-autoclaveable parts		
VN-102	Medicine Cup, non sterile	1
VN-106	Hose System 22 male/15 female	1
VN-109	Filter Membrane	10
VN-111	Luer/Lock Oxygen Hose	1
VN-115	Safety Box with Set of Protecting Plugs	1
VN-116	Luer/Lock Protecting Cap	1
VN-118	Measuring Cup	1
VN-122	Special Face Mask >Children – Size 1 with Expiratory Valve Children/0-1 kg Body Weight	1
VN-123	Special Face Mask >Children – Size 2 with Expiratory Valve Children/1-8 kg Body Weight	1
VN-124	Face Mask >Children – Size 3 with Expiratory Valve Children/8- kg Body Weight	1

Autoclaveable parts - 134°C

VN-101	Filter Shell w/ Valve	1
VN-103	Dome	1
VN-103 komplett	Dome with Sealing Ring, 4 Baffle Plates and Luer/Lock Protecting Cap	1
VN-103 internat.	Dome with Sealing Ring, Blue Baffle Plate and Luer/Lock Protecting Cap	1
VN-104	Exhalation Piece	1
VN-105	Mouthpiece	1
VN-110	Sealing Ring	1
VN-114	Microwave Vapour Sterilizer	1
VN-117B	Baffle Plate Blue - autoclaveable - 3.2 µm	1
VN-117G	Baffle Plate Green - autoclaveable - 2.3 µm	1
VN-117R	Baffle Plate Red - autoclaveable - 3.8 µm	1
VN-117Y	Baffle Plate Yellow - autoclaveable - 4.5 µm	1
VN-117	Set of Baffle Plates (Blue - Green - Red - Yellow)	1
VN-B-109	Tube Extension 22 OD/15 ID, 20 cm long	1
VB 3 Mon.	Autoclaveable Parts for VENTA-NEB®-ir Home Therapy – for 3 Months (incl. 2 Filter Shells, Mouthpiece, Exhalation Piece, Dome with Baffle Pl., Sealing Ring and Luer/Lock Protecting Cap*	1
VB 3 Mon.	Compl. Consumables for VENTA-NEB®-ir Home Therapy – for 3 Months (incl. Autoclaveable Parts, 100 Filter Membranes and 100 Medicine Cups)	1

* The Luer/Lock protecting cap is not suitable for autoclaving or cleaning in the microwave oven

15.0 Warranty

2 years starting from purchase date

16.0 Declaration of Conformity

Manufacturer: NEBU-TEC med. Produkte Eike Kern GmbH

Address: Kreuzfeldring 17
63820 Elsenfeld - GERMANY

Tel.: (+49) (0)6022-610 62-0

Fax: (+49) (0)6022-64 98 12

e-mail: nebu-tec@t-online.de

web: <http://www.nebu-tec.de>

Product Designation: VENIA-NEB®-ir

Model/Type: VN-100/2-2.4 MHz

We herewith declare that the above product complies with the requirements of EC Directive 93/42/EEC.

Applied standards:

Quality System	DIN EN ISO 13485:2003
Electrical Safety	DIN EN 60601-1:1996
Electrical Safety (POMS)	DIN EN 60601-1-4:1996
EMV Standards	DIN EN 60601-1-2:2002
Risk Management	DIN EN 14971:2001

17.0 Warranty Card

See page 27.